

Wish:
Wash, IUSS, Sterilization and High level disinfection insights to improve performance

By
V. Robin O. Novak, RN CIC
Excellentia Advisory Group



Objectives



The participant will be able to:

- ❖ List key steps in the sterilization/ high-level disinfection process
- ❖ Recognize sterilization/ high-level disinfection irregularities
- ❖ Identify cause of instrument discoloration

INTRODUCTION



- Healthy People 2020:
<https://www.healthypeople.gov/2020/topics-objectives/topic/healthcare-associated-infections>
- National action plan to prevent health care-associated infections: Road Map to Elimination :
https://health.gov/hcq/prevent-hai-action-plan.asp?_ga=1.87129296.999411510.1481313804

WASH 

- Most important step in decontamination process
- Removes rather than kills micro-organisms
 - Not microbicidal
 - Remove visible debris
 - Results in low bio-burden; vital to subsequent processes
- Remaining organic residue can inactivate chemical disinfectants or sterilants or protect micro-organisms from destruction
- Debris may become dislodged causing health risks

WASH 

Where do we begin?

- Free of gross soil during surgery/ procedure
- Point of use
- Separate instruments
- Disassemble
- Fully open jointed instruments
- Follow manufacturer IFU

WASH 

- Detergent or detergent/ enzymatic solutions
 - Consult manufacturer IFU
- Blood adheres to surfaces
- Must be mechanically removed
- Proteinaceous blood components become insoluble when exposed to thermal or chemical means

WASH 

Choosing a cleaning agent:

compatible with medical device being cleaned

- nontoxic
- Easy rinse
- Cost effective
- Bio-degradable
- Non-abrasive
- Long shelf life
- Low foaming
- Rapid soil dissolve

WASH 

Method of cleaning depends on device IFU
Manual Cleaning:

- Recommended for delicate / complex devices
- Cleaned under water
- Brush of appropriate size, length, bristle type and material
- Flushed w/ cleaning solution, then flushed with treated water
- Monitor temperature
- Rinse remove debris /detergent residue

WASH 

Mechanical Cleaning:

- Remove soil and micro-organisms
- Mechanical cleaning equipment includes:
 - Utensil washers and cart washers
 - Washer-sanitizers,
 - Pasteurization equipment
 - Washer-disinfectors
 - Washer-decontaminators
 - Washer-sterilizers

WASH



Washer Disinfecter



Removes soil and provides at least 1,000-fold reduction in # of organisms present on instrument by use of hot water (180-200° F) or chemical disinfectant cycle

WASH



Washer Disinfecter –quality check

- Ensure equipment functioning properly
 - Temperature should be monitored and documented.
 - Spray arms should be checked daily
 - Arms free-turning
 - Spray nozzles are not clogged.
 - Strainers cleaned daily or visible debris.

WASH



Mechanical Cleaning:

- Ultrasonic cleaners for fine cleaning
 - Used only after gross soil removed
 - Not disinfection/ sterilization
 - Degassing –follow IFU - document
 - Check IFU to validate U/S cleaner compatible
 - Change cleaning solution before it becomes heavily soiled
 - Thorough rinsing
 - Preventative maintenance per IFU and document

WASH 

Rinsing:

- Removes loosened debris and detergents
- Copious amounts of tap water
- Final rinse:
 - Treated water
 - DO NOT USE SALINE –
 - Salts cause instrument deterioration
 - May interfere with sterilization

WASH 

Verification:

- Visual inspection with lighted magnifier
 - Preventative maintenance per IFU
- Measure organic residue:
 - ATP
 - Residual carbohydrates, glucose, heme

WASH 

Light or Dark Water Spots



WASH 

Light or Dark Water Spots

Cause:

- Air dried instruments
- Mineral content in tap water or condensate in sterilizer

Solution:

- Purified water through reverse osmosis or deionization



WASH 

Rust colored deposits:



WASH 

Rust colored deposits:

Cause:


- Inadequate cleaning or drying
- Inappropriate water used for rinsing
- Incomplete rinse
- Electrolyte deposits underneath chipped chrome plating on stainless steel
- Presence of chlorides

Solution:

- Correct loading of autoclave trays
- Repair and replace the chrome plating




WASH



Blue-gray staining:

WASH



Blue-gray staining:

Cause:

- Chemical solutions
- Cold sterilant solution


Solution:

- Verify type and concentration of chemistries

WASH



Purple-Black Staining



WASH




Purple-Black Staining

Cause:

- Chlorides/detergent in steam lines
- Exposure to saline, blood, potassium
- 2 types of metal placed in u/s together

Solution:

- Clean steam lines
- Verify instrument compatibility



PACKAGING



PACKAGE SELECTION

Sterilization Packaging Characteristics:

- Non-linting
- Barrier to microorganisms
- Resist tearing/ puncture
- Allows sterilant to penetrate, contact item and then be removed
- Suitable for items being sterilized
- Free of toxic ingredients & non-fast dyes
- Complete & secure closure/seal; tamper-evident
- Maintain sterility until pkg opened
- Ease of aseptic presentation
- Cost-effective
- Obtain & keep on file mfr's test data, IFU, care & handling instructions
- P&P based on manufacturer's written IFU

PACKAGING



SELECTION:

- Tensile strength = fabric resistance to breaking when pulled apart
- Higher # = increased resistance to breaking

Wrap Weight	Ave. Grab Tensile Strength
Very heavy weight	110.8
Heavy weight	8.8
Moderate to heavy weight	81.7
Light to moderate weight	65.2
Light weight	45.2
Very light weight	42.1

Source: Cardinal Health

PACKAGING 

PACKAGE PREPARATION

- Clean, dust free and away from fluorescent/ UV light
- Room temperature with humidity between 30-60% for at least 2 hours
- First in, first out rotation
- Examine for defects – discard
- Items being wrapped must be clean/ dry
- Disposable wraps: one and done
- Reusable textile wrap must be laundered between use

PACKAGING 

PACKAGE PREPARATION


- P&P for packaging techniques based on IFU
 - Sequential double-wrapping
 - Simultaneous double-wrapping
 - Wrappers kept snug to prevent low spots
-  Manufacturer IFU for recommended technique

PACKAGING 




© 2015 Mediant

PACKAGING 



ABRASION HOLE
Extreme abrasion with holes in one or both layers of wrap

PACKAGING 



PIN HOLE
A very small hole

PACKAGING 



RUPTURED SEAL


- Long fingernails rupture pkg
- Overstuffed storage can cause pack to become crushed, bent or compressed
- Dropped packages are ALWAYS considered contaminated

PACKAGING 

KNIFE CUT

Cause:
Knife used to open box of wrap


- Straight with clean edges
- Appear intact when pulled together



PACKAGING 

SNAG CUT

Straight cut with a small line of cuts or holes



PACKAGING 


SNAG HOLE

Triangular shaped hole




PACKAGING 

PRESSURE HOLE
A small hole or cut less than 2mm in length



PACKAGING 

PUNCTURE HOLE
Frayed or loose fibers attached to the edges



PACKAGING 

MELTED HOLE
Large round or irregular hole with hard welded edges



PACKAGING 



WET PACK

Related to:


- Steam supply
- Poorly maintained condensation return
- Low steam pressure
- Seasonal changes

PACKAGING 



APPROPRIATE TRAY HANDLING

PACKAGING 



Peel Pack

PACKAGING



Loading

- Stand on edge
- All face same direction
- None touching

PACKAGING



What's wrong with this picture?



PACKAGING



What's wrong with this picture?



Seal is 100%
on plastic

PACKAGING



What's wrong with this picture?



Seal is folded at an angle causing a gap and is 100% on plastic

IUSS



Immediate Use Steam Sterilization:

- Designed for the steam sterilization
- Shortest time between sterilized items removal from the sterilizer and aseptically transfer to sterile field

- Unwrapped instruments
- Abbreviated exposure time & possibly no dry time
- Intended for immediate use; not storage
- NO implantable devices unless emergency
- NOT routine, but exception

• NOT A SUBSTITUTE FOR INSUFFICIENT INSTRUMENT INVENTORY

IUSS



- Packaging system must be compatible with IUSS
- Protect item sterility

Rigid Container systems

- FDA approved
- Suitable for items
- MFR IFU
- Toxic free, low lint, odor free
- Tamper-evident seal
- Allow for ID of contents prior to opening
- Maintain sterility until opened
- Easy to use



IUSS 

- Chemical indicator specific to the sterilization method; used every cycle
 - 1 on inside and 1 on outside
- Used, cleaned, inspected, repaired and maintained per MFR IFU
- Check MFR IFU for limitations
- Packaging must be labeled with instrument, sterilizer, sterilization cycle, personnel and patient
- Staff training/ competency on implementation and yearly



STERILIZATION 

Steam

Definition:

- Complete elimination/ destruction all microbial life
- Physical or chemical processes

Many types: Steam, ETO, Vapor

Advantages :

- Non-toxic ,,
- Cycle easy to control
- Inexpensive ,,
- Rapidly microbicidal ,,
- Least affected by organic/inorganic soils ,,
- Rapid cycle time ,,
- Penetrates packaging, device lumens

STERILIZATION 

Ethylene Oxide (ETO)

Advantages ,,

- Very effective at killing microorganisms ,,
- Penetrates medical packaging
- Cycle easy to control and monitor

Disadvantages ,,

- Lengthy cycle/aeration time
- ETO toxic, carcinogen and flammable
- ETO regulated by some states
- Store cartridges in flammable storage cabinet

STERILIZATION



Hydrogen Peroxide Gas Plasma

Advantages

- Safe and fast
- Leaves no toxic residuals
- Cycle time is 28-52 min
- Heat and moisture sensitive items

Disadvantages

- Cellulose (paper), linens and liquids cannot be processed
- Chamber is small
- Restrictions based on lumen internal diameter
- Requires synthetic packaging
- Hydrogen peroxide may be toxic at levels greater than 1ppm

Sterilization





Does the interior of your facility sterilizer look like this?

Sterilization



Does the interior of your facility sterilizer look like this?



Sterilization 

What's wrong with this picture?



Sterilization 


What's wrong with this picture?



Overloaded, improperly stacked, peel packs lying flat instead of on edge


STERILIZATION 

BOWIE DICK









STERILIZATION 

Biological Indicators



- Gold standard of load sterilization monitoring
- Bacterial spores – *Geobacillus stearothermophilus*
- Requires 2 test with the same lot #
- One in sterilizer. Both in incubator. The one that did not go in the sterilizer should change to a positive result

High Level Disinfection (HLD) 

- Performed on semi-critical equipment
- Requires pre cleaning
- Visual inspection
- Prepared per manufacturer IFU
- Concentration tested
 - Testing strips are tested for validation
- Replaced per manufacturer IFU
- Contact time per manufacturer IFU
 - Appropriate temperature
- Allowed to dry before use
- Stored in a designated clean area to prevent contamination

WISH 

References:

- ANSI/AAMI ST79: AAMI's Landmark Recommended Practice for Hospital Steam Sterilization, 2010
- Rutala WA, Weber DJ. CJD: Recommendations for disinfection and sterilization. Clin Infect Dis 2001;32:1348
- Rutala WA, Weber DJ. Disinfection and sterilization: What clinicians need to know. Clin Infect Dis 2004;39:702
- Rutala WA, Weber DJ, HICPAC. CDC guideline for disinfection and sterilization in healthcare facilities. MMWR. In press.
- CMS infection control surveyor worksheet, exhibit 351

WISH 

QUESTIONS? 

WASH 

Robin O. Novak, RN CIC
 Email: robin@excellentiagroup.com
 Website: www.excellentiagroup.com
 Office: 636-875-5088 ext 101

HEALTHCARE CONSULTING EXCELLENCE!

Excellentia brings healthcare consulting services to client companies within the ambulatory care, ambulatory surgery center and physician practice management niche industries. Excellentia consultants have years of hands-on experience in their respective fields.

Partnering with an industry expert can save you both time and money. Whether you are looking to improve your reimbursement, add a physician to your practice, develop a new surgery center, or need an interim CEO to cover a maternity leave, we have the talent that can help you. Our staff has proven success in their respective field.

Our engagements are flexible in length. We can accommodate projects that are divided into phases as well as long term coverage. Analysis, audits, training or plans of correction are just some of the areas we cover. Our pricing is designed based on the particular needs and can be estimated during our initial conference call once the scope of your project is understood.
