

Sterilization 101: WISH to achieve success

By

Cathy Montgomery, RN, CASC
And

Robin O. Novak, RN, CIC

Virginia
Ambulatory
Surgery
AssociationMaryland
Ambulatory Surgery
Association


Objectives:

The learner will be able to:

- ✓ Identify key steps in the wash, sterilization and high-level disinfection process
- ✓ Recognize sterilization/ high-level disinfection irregularities

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

Sterile Processing Department



STERILIZATION 101: WISH TO ACHIEVE SUCCESS |



STERILIZATION 101: WISH TO ACHIEVE SUCCESS



STERILIZATION 101: WISH TO ACHIEVE SUCCESS



STERILIZATION 101: WISH TO ACHIEVE SUCCESS



WASH

- Detergent/ enzymatic solution
- Mechanical removal
- Visual inspection

STERILIZATION 101: WISH TO ACHIEVE SUCCESS



WASH

Methods of cleaning:

- Manual
- Mechanical
 - Washer disinfector
 - Ultrasonic

STERILIZATION 101: WISH TO ACHIEVE SUCCESS



WASH

RINSING:

- Removes debris/ detergents
- Copious amounts
- Final rinse

STERILIZATION 101: WISH TO ACHIEVE SUCCESS


WASH

Visual inspection:

- Lighted magnifier
- Preventative maintenance
- Validation

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

WASH



STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

PACKAGING

Selection:

- Resistance
- Penetrable
- Suitable
- Tamper evident
- Sterility
- Presentation
- Cost

Preparation:

- Location
- Temperature/ Humidity
- Rotation
- Examination
- Disposable
- Reusable

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

PACKAGING



Techniques:

- Sequential double- wrapping
- Simultaneous double-wrapping




STERILIZATION 101: WISH TO ACHIEVE SUCCESS

PACKAGING




STERILIZATION 101: WISH TO ACHIEVE SUCCESS

PACKAGING



STERILIZATION 101: WISH TO ACHIEVE SUCCESS

PACKAGING



STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

PACKAGING



STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

IUSS

Immediate Use Steam Sterilization:

- Shortest time between removal and transfer to sterile field
- Unwrapped
- Container systems
- Exposure time?
- Implants?



NOT A SUBSTITUTE FOR INSUFFICIENT INSTRUMENT INVENTORY

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

STERILIZATION

Complete elimination/ destruction all microbial life

Steam	EtO	Vapor
<ul style="list-style-type: none"> Advantages <ul style="list-style-type: none"> Non-toxic Easy to control Inexpensive Rapidly microbicidal Disadvantages <ul style="list-style-type: none"> Burns Heat sensitive instruments Wet packages 	<ul style="list-style-type: none"> Advantages <ul style="list-style-type: none"> Heat sensitive equipment Effective Easy to control Penetrates Disadvantages <ul style="list-style-type: none"> Lengthy Toxic, carcinogen, Flammable Regulated Storage 	<ul style="list-style-type: none"> Advantages <ul style="list-style-type: none"> Safe & fast Non-toxic residuals Heat & moisture sensitive Disadvantages <ul style="list-style-type: none"> Paper, linens liquids Small chamber Restrictions Synthetic packaging ? Toxic at certain levels

STERILIZATION 101: WISH TO ACHIEVE SUCCESS

STERILIZATION



STERILIZATION 101: WISH TO ACHIEVE SUCCESS

STERILIZATION



STERILIZATION 101: WISH TO ACHIEVE SUCCESS

HIGH LEVEL DISINFECTION

Inhibits MOST viable organisms except spores and prions

- semi-critical devices
- Visual inspection
- Preparation/ Concentration
- Validation
- Contact time

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

STORAGE & DELIVERY



STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

Quality Assurance Performance Improvement

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |

Variables Affecting the Outcome of Steam Sterilization Process



Ref: Personal Communication, Charles Hancock, President, Charles O. Hancock Associates, Inc.

Four critical variables for steam sterilization to be effective:




Quality Assurance



Monitoring Tools


- Physical monitors
- Equipment monitoring
- Load monitoring
- Pack monitoring
- Exposure monitoring



Physical Monitors



Equipment Monitoring



Load Monitoring - Attest



Pack Monitoring

- Internally monitors pack, trays, peel pouches.
- Verifies the sterilant has penetrated to the point of placement.
- Confirms that specific exposure conditions have been met.



Exposure Monitoring

Quick way to see if a sealed pack has internally been exposed to the sterilization process.



Documentation – 3 years


EXCEPT:
 Alabama
 Florida
 Hawaii
 Maryland
 Mississippi
 Missouri
 Nevada
 N. Carolina
 Ohio
 Virginia



Documentation



Surveyors Corner



Tracking Back



1. Sterilizer load log
2. Cycle number on tape or pouch
3. Call out cycle numbers and record in patients chart

Tracking Back



Transportation Technique



The container must be labeled as biohazard, use a sticker.



Don't Forget To Wear Gloves!!!!!!

REFERENCES

ANSI/AAMI ST79: AAMI's Landmark Recommended Practice for Hospital Steam Sterilization, 2010

- Rutala WA, Weber DJ. CJD: Recommendations for disinfection and sterilization. *Clin Infect Dis* 2001;32:1348
- Rutala WA, Weber DJ. Disinfection and sterilization: What clinicians need to know. *Clin Infect Dis* 2004;39:702
- Rutala WA, Weber DJ, HICPAC. CDC guideline for disinfection and sterilization in healthcare facilities. *MMWR*. In press.
- CMS infection control surveyor worksheet, exhibit 351

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |



Cathy Montgomery, RN CASC
Excellenia Advisory Group
636-875-5088
Cathy@excelleniagroup.com
www.excelleniagroup.com



Robin O. Novak, RN CIC
Excellenia Advisory Group
636-485-2852
Robin@excelleniagroup.com
www.excelleniagroup.com

QUESTIONS?

STERILIZATION 101: WISH TO ACHIEVE SUCCESS |