

healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

33671 Doreka Dr.
Fraser, MI 48026


PRSR STD
U.S. POSTAGE
PAID
MANHATTAN
MAILERS

Endo1501



www.healthmarkgi.com

Visit our new GI website:



Hum!Pak[®] Brushes,
EndoCheck[™] and more

Products for the safe and effective cleaning,
testing, disinfecting and storing of endoscopes,
accessories and supplies.

Endoscopy Products

See Clean Inside a Scope with Endocheck™



Flexible endoscopes are notoriously difficult to clean. More difficult still is verifying that the cleaning has been adequate. Now there is a test for checking the cleanliness of the biopsy channel of the flexible endoscope. The EndoCheck™ is a miniature chemistry kit that is simple to use and interpret. Simply swab the biopsy channel of the scope with the included soft-tipped long probe, clip off the swab into the vial. Mix the activating agent, shake vigorously, wait and then check for a color change. Depending

on the type of test used, a color change indicates that blood residue or protein residue remains in the channel, and should be reprocessed. EndoCheck™ complies with ASTM Guide D7225. Swab diameter sizes from 1.0 - 5.0 mm are available.

Code No.	Description	Price
EDH-110	EndoCheck: Hemoglobin Residue for 1.1mm x 200cm lumen (12 per box)	\$101.55
EDH-200	EndoCheck: Hemoglobin Residue for 2mm x 230cm lumen (12 per box)	\$101.55
EDH-270	EndoCheck: Hemoglobin Residue for 2.7mm x 230cm lumen (12 per box)	\$111.95
EDH-350	EndoCheck: Hemoglobin Residue for 3.5mm x 230cm lumen (12 per box)	\$111.95
EDH-470	EndoCheck: Hemoglobin Residue for 4.7mm x 230cm lumen (12 per box)	\$209.85
EDP-110	EndoCheck: Protein Residue for 1.1mm x 200cm lumen (12 per box)	\$101.55
EDP-200P	EndoCheck: Protein Residue for 2mm x 230cm lumen (12 per box)	\$101.55
EDP-270	EndoCheck: Protein Residue for 2.7mm x 230cm lumen (12 per box)	\$87.05
EDP-350	EndoCheck: Protein Residue for 3.5mm x 230cm lumen (12 per box)	\$111.95
EDP-470	EndoCheck: Protein Residue for 4.7mm x 230cm lumen (12 per box)	\$209.85

Identify Residual Organic Soils with ChannelCheck™




ChannelCheck™ is the first product capable of testing virtually any lumened instrument for residual organic soils, including flexible endoscopes, no matter the channel size. ChannelCheck™ tests for three common organic soils at once: blood, protein and carbohydrates.

Simply flush an individual channel with sterile-DI water and recover that water in a clean (preferably sterile) container. Then dip and swish the test strip in the recovered water for 10 seconds.

Remove the strip, wait 90 seconds and then compare the color of the pads to the results chart. Should any of the pads indicate there is residual soil, reclean the device and then retest.

Code No.	QTY	Price
UCC-101	100	\$103.65


Flexible Endoscope Brushes



Long enough, small enough and flexible enough to clean your flexible endoscopes. Two innovative styles of disposable endoscope cleaning brushes. The Colonoscope/Duendoscope style include the nylon brush safely scrubs the operating & suction channels. The Control Head Brush cleans the opening of the valve bodies.

Code No.	Description	Dia (mm)	Lgth (cm)	Trim (mm)	QTY	Price
00711610	Control Head/Valve Cleaning Brush	1.0	4	10	1	\$6.25
CC-601	Standard Colonoscope/Duodenoscope Color Code: Pink	5.5	230	1.7	12	\$31.20
CC-603	Jumbo Colonoscope/Duodenoscope	9	230	1.7	12	\$31.20
CC-110	Clear Channel Brush Color Code: White	1.2	65		12	\$31.20
CC-172	Tapered Bristle Brush (1.7 - 2.2 mm) Color Code: Purple & Blue	1.7	230		12	\$31.20
CC-374	Tapered Bristle Brush (3.7 - 4.3 mm) Color Code: Orange	3.7	230		12	\$31.20
CC-250	Clear Channel Brush (2.6 - 3.2 mm) Color Code: Green & Yellow	2.6	230		12	\$31.20
CC-601-3770	Standard Colonoscope/Duodenoscope w/Valve Brush Kit	5.5	230	1.7	12	\$54.00
CC-603-3770	Jumbo Colonoscope/Duodenoscope w/Valve Brush Kit	9.0	230	1.7	12	\$54.00
CC-110-3770	Clear Channel Brush w/Valve Brush Kit	1.2	65		12	\$54.00
CC-172-3770	Tapered Bristle Brush w/Valve Brush Kit	1.7	230		12	\$54.00
CC-250-3770	Tapered Bristle Brush w/Valve Brush Kit	2.6	230		12	\$54.00
CC-374-3770	Clear Channel Brush w/Valve Brush Kit	3.7	230		12	\$54.00

Durasponge



The durasponge allows for safe cleaning of surgical instruments and surfaces and it is available in a flat or tubular design.

Code No.	Description	QTY	Price
ESD904101	Tubular - Dry	100	\$62.20
ESD905201	Flat - Dry	200	\$62.20

Endo Parts Bag



The Endo Parts Bag allows you to keep all of the valves for each endoscope together. Detaching the valves before storage keeps the endoscopes from retaining water and therefore prevents contamination.

Code No.	QTY	Price
EBP-46	10	\$93.30

10 Minute Gram Negative Test

Flexible endoscopes are marvelous diagnostic tools. But as we all know, they are a challenge to clean and disinfect/sterilize. Modern machines such as AER's are effective tools for rendering scopes safe for use on the next patient. But what about the endoscope that is not used for multiple days? Is that scope hanging in the closet still safe for use? Well Healthmark has developed a tool kit to help out: the HangTimeKit™. The first step in managing scope storage is to know when they were last reprocessed. For that purpose, the kit includes HTKLabels. Each label is imprinted with the 12 months and with the number 1 - 31, representing the possible 31 days in a month. Simply circle the month and date when the scope is reprocessed and wrap the tag on the scope with the self looping design. When the scope is next used, note the date of last reprocessing in your records. What about a scope that has been stored for more than a couple of days? Is it still safe? For that, the kit includes the HTKTest. This simple and rapid test checks for gram negative bacteria that might have grown out during storage. Utilizing a unique enzyme detection method, the easy to read fluorometer checks for possible colony formation by reading telltale florescence in the recaptured water. If the fluorometer reading is positive for gram negative bacteria, reprocess the endoscope following manufacturer guidelines prior to use.



Clean Cynch Lock



These clean tags can be pulled to varied degrees of tightness. They are ideal when a tight-fitting seal is required.

Code No.	QTY	Price
5030 White Clean	100	\$25.95

Code No.	Description	QTY	Price
403225 HTK	Label with Month and Dates	500	\$150.30
FQ360	Fluorometer (1 each)	1	
HTK-1	Complete starter kit including 1 box of labels, 1 box of gram negative tests & pipettes & 1 each fluorometer	1	
HTK-FV1041	Gram Negative Bacteria Test, 1ml Sterile Pipette (50 per box)	50	\$310.90
403225 HTKKG	Green HangTime Label	100	\$30.05
403225 HTKY	Yellow HangTime Label	100	\$30.05
403225 HTKKB	Blue HangTime Label	100	\$30.05
403225 HTKOK	Orange HangTime Label	100	\$30.05
403225 HTKPK	Pink HangTime Label	100	\$30.05

Transport Trolley



The 2220 Trolley is perfectly suited to transport Healthmark's 2220 round bins, and has the capability to hold up to 5 bins simultaneously. The stainless steel frame is sturdy enough to address the rigors of daily use. Locking casters insure no unwanted movement. This cart and tray system is ideal for transporting clean and dirty scopes to and from reprocessing. By adding Healthmark's green and red tamper evident seals it becomes an excellent process control system for both Clean (green) and Dirty (red) Flexible Endoscopes.

Code No.	Depth	Width	Height	Price
2220 Trolley	32.150	24.350	50.450	\$1000.00

Code No.	Description	Qty	Price
6020 Red	Cynch-Lok	100	\$12.10
6020 Green	Cynch-Lok	100	\$12.10

Round Soaker



Designed by endoscopy nurses, in this 20" diameter tray flexible scopes can be safely coiled for complete immersion in disinfecting solution. The lift-off lid completely covers the tray, limiting fumes. Max temp is 185F.

Code No.	Length	Height	Price
2220	20.000	24.350	\$134.35

Lint-Free Wipes



Manual cleaning of endoscopes and other surgical instruments is not only important, but absolutely necessary. Use these low-linting wipes immediately after use, and prior to automated or manual disinfection. These wipes will help ensure proper reprocessing of instruments.

Code No.	QTY	Price
CC1212	150	\$51.85
CC555	200	\$43.55

Self-Looping Labels



These self-looping labels are sure to be the solution to managing your scope reprocessing and insuring quality patient care. It allows all important information to be written on the label such as department, model number, date time, and pass fail specifications.

Code No.	QTY	Price
404225-351-GI	500	\$233.20

Cool Aids



There are few better ways to beat the heat within the SPD or Operating Room than these lightweight, latex free cooling devices. Simply submerge the head gear or neck band in cold water for 2-3 minutes, and the enjoy hours of cool temps. For additional heat relief pair in one of the Cool Aids vests, which utilize cold packs core body cooling.

Code No.	Description	Price
ECN-001	Evaporative Cooling Neck Band	\$5.60
ECND-001	Evaporative Cooling Neck Band - Wide	\$10.00
ECB-001	Evaporative Cooling Beanie	\$13.75
ECS-001	Evaporative Cooling Skull Cap	\$13.75
ECB-002	USA Evaporative Cooling Beanie	\$13.75
VEST-S	Hybrid Sport Cooling Vest- Small	\$260.00
VEST-M	Hybrid Sport Cooling Vest-Medium	\$260.00
VEST-L	Hybrid Sport Cooling Vest-Large	\$260.00
VEST-XL	Hybrid Sport Cooling Vest-XL	\$260.00
VEST-XXL	Hybrid Sport Cooling Vest-XXL	\$260.00
VEST-XXXL	Hybrid Sport Cooling Vest-XXXL	\$260.00
CPX-001	Replacement Cool Pax Inserts (Set of 4)	\$150.00

PPEs for Endoscopy



The Sleeve Gloves are a unique fusion of an 11 mil nitrile glove and a 4 mil protective sleeve. With an overall length of 28 inches, these gloves provide complete durable hand and arm protection. The sleeve gloves are comfortable, flexible, and lightweight, while the elastic cuff on the sleeve gloves inhibits slippage of the sleeve from the upper arm. These latex-free gloves are a one day use and are resistant to a broad range of hazardous chemicals. Sold in Pairs. **The Drape Long Shield** is ideal for use in departments where maximum personal protection is the goal. This full face shield which provides the best splash protection against exposure to bloodborne pathogens. It adds an 8 inch fluid barrier that extends below the integral full face shield. The barrier fabric tucks inside a sterile gown, providing additional protection to prevent blood from splashing under the shield. **The Decontam Glove** is certified for use in chemotherapy. These gloves are thick enough and tough enough for use in Decontam. This thicker, 15 mil, textured, powder-free latex glove is the decontam glove of choice by professionals who handle potentially hazardous items.



Code No.	Description	QTY	Price
41350	Small Long Sleeve 8ml latex-free glove/4ml sleeve (50 pairs)	50	\$320.00
41450	Medium Long Sleeve 11ml latex-free glove/4ml sleeve (50 pairs)	50	\$320.00
41550	Large Long Sleeve 8ml latex-free glove/4ml sleeve (50 pairs)	50	\$320.00
41650	X-Large Long Sleeve 11ml latex-free glove/4ml sleeve (50 pairs)	50	\$320.00
DSK-1	Small Chemotherapy/Decontam Glove - 15 ml	50	\$16.60
DSK-2	Medium Chemotherapy/Decontam Glove - 15 ml	50	\$16.60
DSK-3	Large Chemotherapy/Decontam Glove - 15 ml	50	\$16.60
DSK-4	X-Large Chemotherapy/Decontam Glove - 15 ml	50	\$16.60
DSK-5	XX-Large Chemotherapy/Decontam Glove - 15 ml	50	\$16.60
18000W-100	Face Shield w/Drape Wide & Long Shield	50	\$112.50
18000-100	Face Shield with Drape	50	\$112.50
16000-100	Face Shield without Drape	50	\$92.50